


The New
VIRACOR
HIGH FREQUENCY VIOLET RAY
GENERATOR

CONTROLLED

ELECTRICITY

for

HEALTH
BEAUTY
VITALITY



MASTER ELECTRIC CO.

SUCCESSORS TO

VIOLET RAY CORPORATION
CHICAGO





Instruction Book

for operating the

VIRACOR

HIGH FREQUENCY VIOLET RAY GENERATOR

IN order to obtain the greatest possible benefit from the use of the VIRACOR, a careful study of the following pages is recommended. In addition to complete instructions, many suggestions for the user are to be found.



TRADE MARK REGISTERED

MASTER ELECTRIC CO.

SUCCESSORS TO

VIOLET RAY CORPORATION
CHICAGO

Printed in U. S. A.



Controlled Electricity . . .

for

THERE is a very well founded theory, that an extremely close relationship exists, between the electrical spark, and the well known "Spark of Life." While the exact relationship is not actually known, numerous experiments and continued laboratory research have proven that electricity, properly controlled and applied to the body by certain means, will give stimulus to the vital organs and help them create a barrier of resistance to the inroads of age, and disease and check loss of vitality.

Controlled Electricity

Controlled electricity is the most modern means of aiding Mother Nature in her wonderful work of renewing health, easing pain and attacking the ravages of disease.

The VIRACOR H.F.1 High Frequency Violet Ray Generator is the final answer to many difficult experiments and long scientific study and research to find the best method of applying controlled electricity for this purpose.

With this wonderful machine, you can aid Nature in your own home. No special experience or knowledge of electricity is necessary. It is simple to operate and easy to handle. You just follow the directions, which come with each machine, and are given in plain non-technical language.

What the Viracor H. F. 1 Does

In the High Frequency Violet Ray the electrical oscillations are of so great a speed that there is only a pleasing warm sensation, no shock or other disagreeable feeling. The Ray is highly invigorating and a great benefit to the entire human body. Through the specially vacuummed glass electrodes, heat and ozone are generated and the electrical energy sent deeply into the body reaching every cell and tissue. The action of these three elements are positive and are a sure method of relieving pain and building up of depleted nerves and weak or sluggish organs.

The viracor H.F.1 High Frequency Violet Ray by its deep penetrating qualities, goes far into the tissues and revitalizes the blood thru cellular massage, thus increasing

circulation, causing a warm flow of revitalized blood to surge thru the body dispelling fatigue, relieving pain and attacking the disease germs which are present in everyone's body. When it is applied to any portion of the body, no pain or any but the most pleasing and soothing sensation is felt.

Not Only For the Sick

The Violet Ray is not to be thought of only as an aid to the sick or injured, but even though you are perfectly well and healthy it will act as a check on reduced vitality and prevent illness and keep you in perfect health.

Use it for the sick or well, men or women, boys or girls, infants or the aged, with the feeling of safety and the assurance that only good can come from its use.

High Frequency Violet Ray treatments and massage are endorsed, prescribed and used by Doctors, Hospitals, Sanitariums, and Beauty and Hair Specialists.



As a Germicide

The VIRACOR H.F.1 has a germicidal action because it generates Ozone. Therefore, frequent application will greatly help to ward off disease and infection and generally tone up the entire system.

For Beauty and Rest

A few minutes a day with the VIRACOR H. F. 1, combined with your favorite facial cream, will keep the skin soft and smooth, and prevent premature wrinkles. The cheeks will have that youthful, healthy glow so much to be desired.



Controlled Electricity . . .

for

If at the end of the day you have that tired and fagged out feeling, use your VIRACOR H.F.1 for a few minutes, and you will be refreshed for an evening of pleasure and contentment with your family, friends, or a book.

Beneficial and Harmless

There is no need to be afraid of this wonderful Ray, it cannot harm you, merely follow the general directions, and the more you use your VIRACOR, the more you will realize the great value and benefits you will derive from its use. Study your own self, and you will find more and more use for your VIRACOR, and soon realize what a wonderful aid to the human race controlled electricity is, when used by means of the VIRACOR H.F.1.

Some of the most important benefits which your VIRACOR will bring you are increased blood circulation, normalized metabolism, soothed nerves, stimulated glandular activity, restored energy, improved digestion, restful sleep, and a general revitalization of the entire system. All this can be accomplished if you will consistently and intelligently use your VIRACOR.



While we do not maintain that the VIRACOR is a "cure-all," we do maintain that the High Frequency Rays which it generates is a wonderful aid in greatly relieving and shortening the duration of pain and disease if used consistently, soon proving, what a boon to humanity it is. The use of High Frequency

Electrical Current is a great benefit to all living matter, it builds up the forces of nutrition and general health, giving that much required "vim and vigor" which is so vital, to be able to cope with the swift progress of life today.

Aids the Work of Nerve Centers

The Spinal Cord is an extension of and is directly connected to the brain. It is encased in the sectional bone structure called the spinal column or the spine.

The spinal cord is the main nerve of the body and from it, branch nerves go to each and every part of the body and its many organs. The nerve network may be compared to an electrical system and just as any trouble on the electrical system is recorded in some manner on the meters of the power houses, any trouble in the body is recorded in the meter of the body, the brain. This recording of trouble is called pain and designates the place where the ailment is located, warning the body that some organ of the body is not functioning properly.

Pain or tenderness at any point along the spine indicates that the organ affected or controlled by that particular nerve branch requires attention. Both the spine and the affected organ should be treated. In most cases, the use of the High Frequency Violet Ray treatments, show exceptional results. The Violet Ray stimulates these nerve centers and the stimulation sent thru the branch nerve track to the affected organ—increasing its activity and causing it to regain its normal tone and resume its duties. When this is accomplished the body is freed from the pain and distress caused by that particular ailment.





About the Treatments

The treatments listed, are just a few tried and tested by physicians and clinics and have proven very beneficial, showing exceptionally good results. After using the Viracor H.F.1 for a while you will develop methods and a technique suitable to your own particular liking.

The treatments need not be of long duration. Three to fifteen minutes treatment at frequent intervals are the most effective.

Under no circumstances should you give any treatment of an orifical nature without first consulting a reputable physician who will give you the proper directions or will treat you himself, and should not be of more than five minutes duration.

Many of the treatments listed are for the aid of the numerous physicians using our products.



Read the "General Instructions" and the "Donts" listed on pages 14 and 15.

While not every case can be entirely cured by the use of the VIRACOR H.F.1 Violet Ray alone, yet practically any ailment can be very materially aided by the use of the VIRACOR Violet Ray as an aid to medical treatment, and many cases can be cured by Violet Ray after medicines have failed.

Suggested Treatments

ALOPECIA or FALLING HAIR—Use Electrode No. 3, mild current at first, increasing

with treatments. Pass the comb over the scalp for five minutes each day. It is suggested that a hot oil treatment be used in connection. The same treatment for baldness and restoration of grey hair to natural color.

ACNE (PIMPLES and BLACKHEADS)—Use Electrode No. 1. Treat entire surface of affected parts for about five minutes. If skin is damp either apply talcum powder or treat through a handkerchief. During period of treatments watch your diet and avoid rich and greasy food. Take a mild laxative.

ABSCESS—Use Electrode No. 1, keep in close contact with skin, using enough current to get full benefit of heat and ozone generated. Treatment three to five minutes once or twice a day.

ANEMIA—Use Electrode No. 1. Treatment eight to ten minutes daily Treating over entire spine, chest, abdomen and back. Use the Ozone Generator No. 28 in connection.

ASTHMA—Apply Surface Electrode or Condenser No. 2 over the chest and throat, treating for about five minutes, following with treatment of Ozone Generator Electrode No. 28.

BACKACHE — Use Electrode No. 1 over sore part of the back for about five minutes, using rather a strong current. (Backache may be caused by weakness or disease of some internal organ, and it is suggested that a Physician be seen).



BARBER'S ITCH — After cleansing the infected area with a mild antiseptic wash, treat with surface Electrode No. 1 in close contact with skin, followed by short sparks for about a minute. Continued treatment is usually required as this condition is a most stubborn one.

BLADDER TROUBLE (Cystitis)—With the use of Electrode No. 25 and urethral Electrode No. 21 remarkable results have been obtained. Treatment five minutes daily, using a good antiseptic lubricant on the electrode.

BRONCHITIS—Treat the same as for Asthma.

BRAIN FOG—Use Electrode No. 1 over the forehead, eyes and back of neck and head with a spark just strong enough to cause a good tonic reaction, followed by two or three minute treatment with Ozone Generator No. 28.

BREAST DEVELOPMENT (Flabby or drooping Breasts)—Use Electrode No. 1 from five to eight minutes daily. Cover breast with cocoa-butter, then move your VIRACOR lightly over each side from the neck downward and also from armpit to armpit under each breast, pressing lightly upward. Treat 10 minutes daily. Remarkable results may be expected.

BRUISES—Use Electrode No. 1 over affected area, five to ten minutes starting with medium and ending with a strong current. Use an ointment in connection.

BOILS—Treat the same as for acne.

BURNS—Apply Electrode No. 1 in close contact with burned area, using a mild current for a short period, two or three time daily. Use a good ointment in connection.

CARBUNCLE—Use Electrode No. 5 in loose contact with inflamed area about five to ten minutes daily.



Controlled Electricity . . .

for

CATARRH (Nasal)—Use Electrode No. 1 over external area for five minutes on each side, followed by treatment with Ozone Generator No. 28.

CATARRH OF WOMB—Use vaginal Electrode No. 24 or 26, five minutes daily. Use a good lubricant. Mild current.

CHILBLAINS—Use Electrode No. 4. Wash feet thoroughly in warm water, thoroughly dry, apply carbolated vaseline, treatment about ten minutes on each foot.

COLD HANDS or FEET—Use Electrode No. 2 or 4, from five to ten minutes two or three times daily. Use strong current.

COLD IN HEAD—Same treatment as for nasal catarrh.

COLD IN LUNGS—Same treatment as for asthma.

CONSTIPATION—Use Electrode No. 1, massage the entire abdominal region at first directly in contact with the skin, later followed by a short spark over the same area. Ten minutes daily at first, ten two or three times a week. In particularly stubborn cases the use of rectal Electrode No. 17 properly lubricated, three to five minutes.

CORNS—Use Fulguration Electrode No. 25 for 10 to 25 seconds daily. Care should be used to apply spark on the hard tissue only and not on the surrounding area, then Electrode No. 18 for three minutes.

DANDRUFF—Hot oil shampoo before starting treatments. Use Electrode No. 3, the same as a comb, treating three to five minutes every two or three days (do not use in connection with hair tonic containing alcohol).

DEAFNESS, EARACHE, ETC.—Apply ear Electrode No. 11 about three minutes daily, followed by use of electrode No. 1 in back of the Ear and over side of face.

DIABETES—Use Electrode No. 1 or 2 over the abdomen and spine through a thin cloth or towel, apply strong current about ten minutes daily. Complete results should not be expected for several months. Diet must be observed, and it is suggested that you consult a competent physician.

DYSPEPSIA—Use Electrode No. 1 or 2 over the stomach and entire abdominal region, using a strong current for about ten minutes daily.



ECZEMA—Use Electrode No. 1 or No. 4. A handkerchief laid over the surface will cause a spark and is preferable. If the skin itches hold the electrode slightly away from the skin. This will both relieve and gratify. Excellent results have been obtained by giving the patient as much current as can be endured comfortably.

EYE DISEASES—Iritis, retinitis, strophy of optic nerve, conjunctivitis, trachoma, glaucoma, incipient cataract, paralysis of ocular muscles, intra-ocular hemorrhage have all been greatly aided by high frequency currents. Close eyelids and apply a mild current. Use Electrode No. 5 or No. 7 for short periods, not exceeding three minutes every other day.

ENTERITIS—Give tonic applications with Electrode No. 1 over the area of the abdomen. Fifteen minutes of treatments daily. In chronic cases, excellent results can be obtained by giving treatment twice daily with the condenser Electrode No. 2.

**EROSIONS OF THE CERVIX
UTERI**—Use Vaginal Electrode No. 24. Excellent results have been obtained. Bring tube in contact with the Cervix, giving ten minutes daily, using antiseptic lubricant on the Electrode.

FELONS—Treat with No. 5 Electrode and a mild current in close contact daily from four to six minutes.

FEMALE TROUBLES — Use Electrode Nos. 15, 24 or 26. Using an antiseptic lubricant before insertion. Turn the current on after insertion of electrode and turn the current off before removing to avoid unnecessary sparking. A medium current for five minutes is applied. Repeat treatment as circumstances require. During menstruation an excessive flow of blood may result, as high frequency currents draw the blood.

FRECKLES—Cover surface with gauze and treat daily with a medium current of from five to ten minutes' duration. Results are quite apt to be very slow and patience is required. Use a good facial cream in connection.

FROST BITES—Treat affected part with No. 1 Electrode, using instructions outlined for the treatment of chilblains may also be adapted to the treatment of this ailment.



GOITRE—The absorption of tissue being necessary persistence is required. Treat with external throat Electrode No. 4 in close contact with the skin and follow with moderately sharp spark over the surface and sides of the protuberance. Ten minutes daily.

GOUT—Use Electrode No. 1 or No. 4. This condition may be helped by a continued massage over the area with a mild current. If pain is felt at first begin with short treatment and increase. Apply spinal Electrode No. 6 to the spine to stimulate the system. Then gently massage with Electrode No. 1 or No. 4 to the painful region with a medium current, keeping the electrode in contact with the skin. Pain may be noticed at first, but results are certain if treatment is continued. Do not treat during acute attack.

GRIPPE (Influenza) — Treat spine, solar plexus, over eyes and sides of nose, using Electrode No. 1, No. 5 and No. 6. Then breathing through Ozone Generator No. 28 gives almost immediate relief. Follow with treatments over side of nose, and in some cases with Nasal Electrode No. 19 internally. Treatment five to ten minutes.

GONORRHEA — GLEET — In sub-acute or chronic forms, marked beneficial results have been obtained by the use of the urethral No. 21 or vaginal No. 24 Electrodes. Apply any antiseptic lubricant on tubes before inserting. Give seven minute treatments daily.

HAY FEVER—Apply Electrode No. 1 or No. 5 over the nose and No. 6 on the spine. Treatment must be of short duration but frequently repeated. Use the Ozone Generator No. 28 freely, it being advisable to Ozonate the entire room with it. For best results treatments should begin two or three weeks before period of attack.

HEADACHES—Use Electrode No. 1 and No. 5 or No. 7. Headaches are from varying causes. If relief is not had from applying to the seat of pain, treat the spine and stomach area. Keep the Electrode in loose contact and prolong the treatment until relief takes place, which usually takes from five to ten minutes.



HIVES AND RASH—Use Electrode No. 1 with medium current strength over the affected area. Use a good ointment before applying the electrode; three to five minutes. Treat daily.

HEMORRHOIDS (Piles) — Use Electrode No. 23 with mild current for about five minutes, care being observed that the electrode is properly lubricated and that current is not turned on until after insertion.

HERPES—Zoster—Use Condenser Electrode No. 2 in the erupted area. Three to five treatments are usually sufficient. Ten minutes daily.

INFLUENZA—Use Electrode No. 1 and No. 6. Treat the spine and solar plexus. Also over the eyes and the sides of the nose. Treatments to the spine are made through a towel for five minutes. Also use Ozone Generator No. 28 freely.

INSOMNIA—Use Electrodes No. 1 and No. 4. Cover back of head, neck, spine and over the eyebrows. Use a strong current five minutes daily before retiring.

LARYNGITIS—Use Electrode No. 4. Either loose contact or mild spark will greatly aid in treatment of this nature. Ozone inhalations are advisable and most beneficial.

LEUCORRHEA—Follow the directions given under female troubles. Also use Electrode No. 6, moving it up and down the spine slowly. A strong current must be used for about five to six minutes daily. Then use vaginal Electrodes No. 24 or 26, lubricating well before inserting. Use a medium current for five minutes and repeat as circumstances require. Current should be turned on after insertion of electrode. Use antiseptic douches in connection.

LUMBAGO—Use Electrodes No. 1 or No. 4. Treat over painful area until all pain is dispelled. Use as strong a current as can be tolerated. Apply the current through the clothing with as strong a current as can be tolerated. Repeat treatment twice daily, until pain is much relieved. Keep the electrode moving and give long treatment about 15 minutes.

MUMPS—Use Electrode No. 1. Treat the swollen parts with a medium current, keeping it in contact for about five to eight minutes. Repeat as necessary, once or twice daily.

NERVOUSNESS—Use Electrodes No. 1, No. 4 or No. 6. Treat the spine, back of head and neck with a mild current. Immediate relief may be expected. Use a weak current at first and increase the strength gradually; five to eight minutes of treatment once or twice daily.

NEURALGIA—Use Electrode No. 1. Apply to seat of pain, after applying lubricant, raising electrode occasionally to produce spark. Pass back and forth in light contact with the skin, giving ten minute treatments.

NEURITIS—Use Electrode No. 1 or No. 4. The first few treatments usually increase the pain, after which relief is felt, and the current may be increased and a short spark given. Several treatments a day are advisable. Use light current at first, and gradually increase same. Apply locally or to entire body, depending upon the condition. Five minute treatments.



Controlled Electricity . . .

OBESITY—Use Electrode No. 1 or No. 2. Apply to affected parts consistently. Begin with a medium current and increase to full power. Rigorous diet, exercise and personal hygiene must, of course, be observed. Fats will be consumed by increase of oxidation. Treatment 15 minutes twice daily.

PARALYSIS—Use Electrodes No. 1 or No. 2. Treat paralyzed muscles, employing a current strong enough to produce a half or three-quarter inch spark. Part of the time keep the applicator in contact with the skin and part of the time raise it above the surface to get effective spark. Treat as directed by physician.

PROSTATIC DISEASES—Use Electrode No. 22. Orifically and electrode No. 1 or 6 on spinal area. This most distressing ailment, common to so many men of middle age, is one in which the high frequency current is almost a specific.

Medical authorities agree that more than 60% of all men past middle age (many much younger) are afflicted with a disorder of the prostate gland. Aches in feet, legs and back, sciatic pains, and impotence are some of the signs. No man should approach the prime of life content to regard these pains and conditions as inevitable signs of approaching age. The VIRACOR may be used to restore youthful health and vigor and to restore the prostate gland to its proper functioning.

In treating the Prostate through the rectum place the patient on one side in the Sims' position with the knees well drawn up. Electrode No. 22 well lubricated, is then inserted in the rectum about six inches with the depression in the tube turned toward the front or anterior wall of the rectum so the gland rests into depression for best results. In Prostatic diseases of all kinds the High Frequency current has proven most efficacious. The duration of each treatment is three to five minutes. Extraordinary results are accomplished in all of Prostatic diseases. Enlarged Prostate should show improvement from the very first—about two to three treatments a week. Treatment by or under direction of a physician is advised. Same treatment for Senile hypertrophy, impotence, etc. Two or three treatments a week then skip a month.

PLEURISY—Use condenser Electrode No. 2 over the front and back of the chest until a marked reddening of the skin indicates a counter irritation. Two or three daily treatments are necessary in acute cases and one daily in chronic cases. Duration 5 to 10 minutes.

RHEUMATISM—Use Electrode No. 1, No. 2 or No. 4. High Frequency currents are of exceptional value in muscular and in chronic articular rheumatism. The pain will be relieved by application of High Frequency massage, using lubricant, if necessary. Follow with fairly strong spark until pronounced redness is produced. This treatment persistently followed has in many cases resulted in a cure. While some relief may be obtained in one treatment, but permanent cure requires much time and patience. Treatment 10 to 15 minutes daily.

RINGWORM—Use Electrode No. 1 and apply short sparks. A medium current is sufficient. Move the electrode so as to cover the entire affected area. Treat two or three times a day, then skip a day.

SCARS—If the tissue is hardened or callous its removal may be facilitated by the application of a mild spark through Electrode No. 25, three to five seconds. Then use Electrode No. 1, by which a softening effect can be produced, treating same five to seven minutes daily. Many treatments are usually required.

SKIN DISEASES—Use Electrode No. 1. Treat over the affected tissues with electrode held in close contact for from five to fifteen minutes. Follow by moderate spark three to five minutes. Treat once daily. If the skin is damp use talcum powder so that electrode will not stick. See treatment for Acne and Eczema.

SORE FEET AND STONE BRUISES—Use Electrode No. 4. Applications in close contact with the skin brings quick relief. Use lubricant and treat for five to ten minutes daily.

SORE THROAT (Tonsilitis)—Use Electrode No. 4. Hold directly against the skin with fairly strong current. The use of the Ozone Generator No. 28 with mouth applicator for several minutes following will be of great benefit.

SPRAINS—Use Electrode No. 1 or No. 4. The pain and soreness following a sprain may be alleviated by treatment first with the electrode in contact with the skin and later by a mild spark, gradually increasing in strength until a hyperemia is produced.

STIFFNESS OF JOINTS AND MUSCLES—Use Electrode No. 1 or No. 4. This condition is generally due to the over-employment of certain muscles. The resulting stiffness and soreness may be relieved, using the same method as explained in Neuritis.

TUBERCULOSIS—In the incipient type of tuberculosis, exceedingly good results can be obtained and tubercule bacilli reduced by the use of the Ozone generator No. 28, filling same with one-half dram of oil of eucalyptol, and one dram of oil of pine, ten minims of beechwood creosote, using the same morning and night. Use Electrode No. 1 over the area of the chest, this will also be found beneficial, using same ten minutes daily.

WARTS AND MOLES—Use Electrode No. 25. These disfigurements, no matter of how long standing, may be successfully treated by the application of a spark. The greatest care must be exercised that these sparks do not come in contact with the surrounding tissues. Permanent results may be expected, and the most disfiguring warts or moles removed by these means.

WHOOPING COUGH—Treat same as Asthma.

STIFF NECK—TORTICOLLIS—Use Electrode No. 1 or No. 2 over the area of the neck and down the spine daily for ten minutes, relief can be effected.


WRINKLES—Are commonly caused by using a given set of facial muscles more than normally. They can be removed by applying Electrode No. 1 or No. 5 with a rotary, massaging movement.



Controlled Electricity . . .

fo

The VIRACOR H. F. 1 Generator



Made of the finest materials obtainable scientifically designed, expertly manufactured and carefully assembled.

It is pleasing to the eye, easy to handle, safe and simple to operate.

EQUIPMENT:

8 feet of Art Silk Cord and Bakelite plug, attached, ready to plug into your light socket or base board receptacle. 3 glass treatment electrodes.

No. 1—General utility electrode, for use on head, face, or body.

No. 3—Comb type electrode, for use on scalp or spine.

No. 4—Throat electrode, for use on throat, arms, limbs and other curved parts of the body.

Comes packed in a beautiful Lacquered Leatherette covered case, lined with royal purple Art Satin. Has neat compartments for storing the electrodes and instruments.

Price Complete — \$12.50

Equipped for standard current (105 to 120 volts A.C. or D.C.) (Can be supplied for other voltages at a small advance in price.)

The cost of operating the Viracor H.F.1 Violet Ray Generator is negligible as the electric current consumed amounts to less than one cent per hour.

General Instructions as to Use and Care of Instrument

1. See that Knob is turned to OFF position. (Turn gently to left as far as possible.)
2. Insert Glass Electrode in holder at small end of the instruments, and see that the Electrode is approximately body temperature.
3. Place Plug in receptacle or socket and turn on current.
4. Turn knob to right until a mild current is flowing through the Glass Electrode.
5. Apply and withdraw Electrode with a quick movement direct to part of body to be treated. (This will prevent the patient from feeling the discharge which occurs when the Electrode is brought in contact with or withdrawn from the body.)
6. Increase flow of the High Frequency Ray to proper strength by turning knob to the right.
7. Follow instructions as given in the book.
8. All Electrodes for internal (orificial) treatments should be thoroughly lubricated with vaseline or other lubricant, and after insertion should be kept constantly in motion to prevent the glass electrode from sticking to the membranes. Turn current on after inserting of Electrode.
9. Under no circumstance should internal (orificial) treatments be taken or given without first consulting a physician and have him either give the treatment or demonstrate the exact course you should follow. This also applies to any serious disease, infection or any acute condition. (A number of the treatments described are given for the use of the physician only.)
10. When treatment is completed turn the knob to OFF position (to the left), remove plug from socket or receptacle, clean and sterilize electrode and return electrode and instrument to the case.
11. It is not the purpose of the VIRACOR High Frequency Violet Ray Generator to supplant the services and attention of a competent physician, but is to be used in conjunction with and in many cases under his supervision and direction.
12. If treatment is of more than 15 minutes' duration, treat for 15 minutes and then allow the generator to cool off for 10 minutes before continuing the treatment.

DONT'S

1. DON'T use instrument near a radiator in the bathroom, with wet hands, or in connection with a liniment, tonic or other medicant containing alcohol, or other highly inflammable ingredients.
2. DON'T make contact with the glass electrode on a metal surface while the current is turned on.
3. DON'T twist or kink the cord.
4. DON'T wait until the cord is entirely worn out, have it repaired when it shows signs of wear.
5. DON'T attempt to repair the instrument yourself. Should it become out of order, send it back to the factory prepaid. A small charge will be made for repairs and adjustments *after* the guarantee has expired.
6. DON'T fail to have glass electrode well lubricated when giving internal (orificial) treatments.
7. DON'T try to be your own physician. If you have an acute or exceptionally infectious disease or ailment, see your physician and he will prescribe for you, and very likely advise you as to how to use your VIRACOR to the best advantage in your particular case.

Ozone Generator



This Electrode is designed for direct medication and treatment of nasal and throat ailments. It generates large quantities of ozone. A very faint blueish vapor, composed mainly of oxygen, a gas in which few germs can live. This gas is formed by passing high frequency electricity thru the air and in its self, is very beneficial. Together with the inhalant exceptional results may be obtained in cases of Catarrh—Asthma—Hay Fever—Bronchitis—Influenza and in fact almost any condition affecting the respiratory organs. The vapor reaches every cell destroying germs and soothing the inflamed tissues.

Attach glass Nasal (D) or Mouth Applicator (E) to rubber hose (B) then attach other end of rubber hose to part designated as No. 1A then attach rubber bulb (C) to part designated as No. 2A, then remove the cork No. 3A, and using the glass funnel (F) fill the chamber (1/5) one-fifth full of inhalant (G) replace the cork and then insert the generator (A) into your VIRACOR H.F.1 Violet Ray machine. Then the

current is turned on. The left hand is then put around the H.F.1 Generator while the other hand actuates the pressure bulb at convenience. The nostrils are placed over or right against the nose-piece and the ozone is forced down through the nose into the lungs by pressing the bulb. The inhaling of the ozone should be deep and the air completely exhaled before taking another breath. Inhalation may also be made through the mouth by changing the nose-piece for the mouth-piece.


The inhalant solution contained in the upper chamber of the generator washes the ozone as it bubbles through it, allowing particles of it to be carried along in the form of vapor. This inhalant solution is a splendid combination of oils which have a soothing effect upon the membranes and greatly assists in healing.

No. 28—Ozone Generator comes equipped with the Generator Electrode, Glass Nasal attachment; glass mouth attachment; glass funnel; rubber connecting hose; rubber pressure bulb and a one ounce bottle of specially prepared inhalant.

COMPLETE

\$12.50

VIRACOR ELECTRODES



The VIRACOR Glass Electrodes are made from the finest, specially made glass and are specially annealed and vacuumed—scientifically shaped by expert glass blowers.

- No. 1—General Surface type Electrode.....\$1.00
For facial massage and general body treatment. Is supplied with all VIRACOR High Frequency Generators.
- No. 2—Condenser Electrode...\$3.00
Gives an exceptionally strong current. Used for deep seated ailments such as rheumatism or lumbago.
- No. 3—Comb or Rake Electrode.....\$2.00
For scalp and hair treatments and may be used on the spine.
- No. 4—External Throat Electrode.....\$1.00
For throat or other curved parts of the body.
- No. 5—Single Eye Electrode...\$1.25
- No. 6—Spinal Electrode.....\$1.50
For all spinal treatments.
- No. 7—Double Eye Electrode...\$2.00
- No. 8—Small Cataphoric Electrode.....\$2.25
Has small indentation to hold cotton and liniment.
- No. 9—Large Cataphoric Electrode.....\$4.00
Has place for cotton and liniment for lumbago, rheumatism, sprains and wherever a liniment is used.
- No. 10—Spatula or Tongue Electrode.....\$2.00
- No. 11—Insulated Ear Electrode.....\$2.25
For ear troubles of all kinds.
- No. 12—Special Metal Electrode.....\$1.25
Delivers a very strong current and used for indirect treatments.
- No. 14—Urethral Electrode...\$1.00
- No. 15—Plain Vaginal Electrode.....\$1.00
- No. 16—Plain Internal Throat Electrode.....\$1.00
- No. 17—Plain Rectal Electrode.....\$1.00
- No. 18—Corn Electrode.....\$1.25
- No. 19—Insulated Nasal Electrode.....\$2.00
- No. 20—Insulated Internal Throat Electrode.....\$2.00
- No. 21—Insulated Urethral Electrode.....\$2.00
- No. 22—Insulated Prostatic Electrode.....\$2.00
- No. 23—Insulated Rectal Electrode.....\$2.00
- No. 24—Insulated Vaginal Electrode.....\$3.00
- No. 25—Cantery Electrode...\$2.00
For removal of warts, moles and other growths.
- No. 26—Special Perforated and Insulated Vaginal Electrode.....\$3.00
Used for ovaritis, leucorrhea, vaginitis and all inflamed conditions of the vagina. Also as a mild germicide. The perforations allow the ozone to be generated and liberated.

Guarantee

ALL VIRACOR PRODUCTS
are guaranteed for a period of
one year against electrical, me-
chanical, and material defects,
providing ordinary care is taken
in their use and handling (Glass
electrodes, electric cord and plugs
excluded).

VIOLET RAY CORPORATION

CHICAGO

The Master  Electric Co.

2309 Archer Avenue
CHICAGO, ILL.

